


# Tubes

## Standard Tubes



- Maximum service temperature of PFA is 500°F/260°C
- Tubes can be autoclaved
- PFA tubes are non-breakable
- Closures sold separately

	Volume	Part Number	Top Opening	Interior	Opening Diameter x Height (mm) (with closure)	Accepts Closure	
						Part Number	Description
A	5 ml	210-005-21	Open	U	11.9 x 74.9	NA	-
B	14 ml	210-014-21	Open	U	16.0 x 99.8	NA	-
C	15 ml	210-015-20	Threaded	U	16.0 x 113.0	600-018-01	Plain threaded
D	15 ml	210-015-30	Threaded	V	16.0 x 129.3	600-018-01	Plain threaded
E	26 ml	210-026-20	Threaded	U	21.6 x 116.1	600-024-01 600-024-23 600-024-71	Plain threaded (2) 1/8" ports Recessed
F	32 ml	210-032-20	Threaded	U	20.1 x 153.4	600-024-01 600-024-23 600-024-71	Plain threaded (2) 1/8" ports Recessed
G	50 ml	210-050-20	Threaded		28.4 x 122.7	600-029-01	Plain threaded
H	50 ml	210-050-30	Threaded	V	28.4 x 119.4	600-029-01	Plain threaded
I	60 ml	210-060-20	Threaded	U	25.9 x 168.7	600-028-01	Plain threaded



Recessed closure with Nomex® tab (tabs sold separately)

V = Conical Interior    U = Rounded Interior



Dimensions available online for download



## Block Digestion Tubes



- PFA tubes compatible with block digestion systems
- Designed for ultratrace applications when traditional polypropylene digestions tubes are not clean enough
- Maximum service temperature of PFA is 500°F/260°C
- Thin walled for maximum heat transfer
- PFA closures sold separately
- Can also be used for long term sample storage

Volume	Part Number	Top Opening	Interior	Base Diameter x Height (mm) (with closure)	Accepts Closure	
					<b>Part Number</b>	<b>Description</b>
<b>A</b> 4 ml	210-004-20	Threaded	U	15.2 x 46.2	600-023-01	Plain threaded
<b>B</b> 10 ml	210-010-30	Threaded	V	26.2 x 40.6	600-033-01	Plain threaded
					600-033-23	(2) 1/8" push-in side ports
					600-033-26	(2) 1/8" ports, one top, one side
					600-033-71	Recessed
17 ml	210-017-20	Threaded	U	26.2 x 40.6	600-033-01	Plain threaded
					600-033-23	(2) 1/8" push-in side ports
					600-033-26	(2) 1/8" ports, one top, one side
					600-033-71	Recessed
<b>C</b> 25 ml	210-025-30	Threaded	V	26.2 x 78.5	600-033-01	Plain threaded
					600-033-23	(2) 1/8" push-in side ports
					600-033-26	(2) 1/8" ports, one top, one side
					600-033-71	Recessed
<b>D</b> 55 ml	210-055-11	Open	—	30.2 x 109.5	-	-
<b>E</b> 60 ml	210-060-70	Threaded	—	29.5 x 132.6	600-033-01	Plain threaded
					600-033-23	(2) 1/8" push-in side ports
					600-033-26	(2) 1/8" ports, one top, one side
					600-033-71	Recessed
<b>F</b> 75 ml	210-075-11	Open	—	34.5 x 104.4	600-038-89	Push-in plug
<b>G</b> 150 ml	210-150-11	Open	—	38.1 x 218.4	-	-



V = Conical Interior    U = Rounded Interior    — = Flat Interior

## Microcentrifuge Tubes

See Microcentrifuge Vials, page 44.



Dimensions available online for download